
MILLIPORE

MILLI-Q 185 PLUS

User manual

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MILLIPORE

Copyright

Copyright © 1990, all rights reserved, Millipore Corporation.
PF04517, Rev. A , 10/93.

Trademarks

Millipore is a registered trademark of Millipore Corporation.
MILLI-RO, MILLI-RX, DURAPORE are trademarks of Millipore Corporation
ANATEL A 100 P is a trademark of ANATEL INSTRUMENT Corporation.

• Table of contents

Table of contents	i	Chapter 5 : Troubleshooting	11
		Troubleshooting table	11
Chapter 1 : Introduction	1	Chapter 6 : Appendix	15
General description	1	Appendix 1 : Main fuse replacement	15
Description of function	1	Appendix 2 : Main P.C. board, fuse replacement	15
Schematic	1	Appendix 3 : Replacement of fuse on U.V. Lamp P.C. board	15
Specifications	2		
Chapter 2 : Installation	3	Chapter 7 : Ordering information	17
Unpacking	3		
Positioning of MILLI-Q 185 PLUS	3	At back of document	
Installation	3	Chapter 8 : Figures	
Inlet water connection	3		
Installation of purification pack Q-PAK	4	Chapter 9 : Additionnal information	
		Technical assistance	
Chapter 3 : Operation	5	Warranty	
Standard modes of operation	5		
Additional displays	6		
Use of keyboard	6		
Electrical connections	7		
Startup	7		
Installation of final filter-MILLIPAK 40	7		
Withrawing water	8		
Chapter 4 : Maintenance	9		
Maintenance schedule	9		
Replacement of Q-PAK purification element	9		
Replacement of final filter MILLIPAK 40	9		
Replacement of the U.V. lamp	9		

1 INTRODUCTION

.....

• General description

The MILLI-Q 185 PLUS has been developed specifically to produce Ultrapure water with the lowest level of organics currently attainable. This water quality can be produced using a laboratory grade feed water prepurified by Reverse Osmosis, or distillation to feed the system.

The system is designed to produce up to 1.5 litres/minute of Ultrapure water on demand.

This water quality should be used for all critical applications requiring the lowest levels of organics : e.g., Critical H.P.L.C., G.C./M.S., Langmuir Blodgett monolayers, I.C.

The MILLI-Q 185 PLUS System main components are shown on figure 1 as follows :

- (A) Water purification module - Q PAK.
- (B) Protected point of use.
- (C) Point of use, final filter, 0.22 µm membrane. MILLIPAK 40 unit.
- (D) Point of use 3 - way valve.
- (E) Control panel.
- (F) Q-PAK locking handle.
- (G) Indicator lights.
- (H) Liquid Crystal Display.
- (I) Keyboard.
- (J) Feed water inlet connection.
- (K) Removeable plug - blocking the socket for the wall mounting bracket arm.
- (L) Knurled screw locking rear cabinet in place.
- (M) Female Jack socket for connection of remote control device.
- (N) Power inlet (mains cord).
- (O) Main fuse holder.
- (P) ON / OFF switch.
- (Q) Heat sink.

• Description of function

The prepurified feedwater (Reverse Osmosis, or Distilled) is driven by a built in pump, through a chamber where the feedwater is exposed to ultraviolet light (U.V.) of low wavelength (185 nm + 254 nm). This initiates oxidation of dissolved organic compounds remaining from the prepurification process, (Reverse Osmosis or Distillation). The U.V. irradiation also displays a bacteriocidal action.

The U.V. treated feedwater is then passed through the Q-PAK purification element which removes the ionic compounds plus the by products of U.V. photooxidation.

A Built-in three way valve allows recirculation of the water until the required water quality is achieved. When the required water quality is reached, turning the three way valve allows water to pass through the MILLIPAK 40 final filter unit at the point of use.

• Schematic

The MILLI-Q 185 PLUS system operating schematic is shown on figure 2 as follows :

- (A) Inlet solenoid valve.
- (B) Pump.
- (C) U.V. lamp in external housing.
- (D) Water purification pack - Q-PAK.
- (E) Resistivity measurement cell to measure product water quality.
- (F) Three way valve.
- (G) Final filter Millipak 40 (0.22 µm).
- (H) Check valve.

• Specifications

Dimensions	
Height	495 mm (19.5 in)
Width	297 mm (11.7 in)
Depth	433 mm (17 in)
Operating weight	17 kg
Electrical specifications	
230 V/50Hz - 180 VA	Fuse 1 Amp slow blow
Hydraulic specifications	
Feed water tubing	External diameter 8 mm Maximum length 3 M
Feed water pressure	Minimum - Positive head (0.3M) 0.03 bar (0.5 psi) Maximum - 1.0 bar (15 psi)
Feed water flowrate	1.5 Litres / min
Performance	
Product water flow rate	≤ 1.5 litres / min
T.O.C. (Total Oxidisable Carbon)	< 5 ppb typical value by ANATEL A 100 P

2 INSTALLATION

.....

• Unpacking

On unpacking the MILLI-Q 185 PLUS you should find the following items, as shown on figure 3.

- (A) MILLI-Q 185 PLUS unit.
- (B) Power cord .
- (C) Point of use mounting cover.
- (D) Users manual.
- (E) Adaptor, 1/2" Female - 8 mm (0,3 in) tubing with built in screen filter.
- (F) Teflon tape.

Following items to be ordered separately

- (G) Water purification Pack - Q-PAK 1.
 - (H) Final filter unit - Millipak 40 (0.22 µm).
- Items G + H = Cat N° CPMQK05R1.

• Positioning of MILLI-Q 185 PLUS

The MILLI-Q 185 PLUS unit may be placed on a laboratory bench, or wall mounted.

To mount the unit on the wall, order the wall mounting bracket (Cat N° ZFMQ00WMB) and verify that the wall is able to support the operating weight of the unit (17 Kg).

Figure 4 shows the connection of the MILLI-Q 185 PLUS with different feedwater delivery options (reservoir or pressurized line).

NOTE : The isolating valve (on the pressurized line), the pressure regulator and pressure gauge shown on figure 4 are not supplied with the system.

• Installation

INLET WATER CONNECTION (FIGURE 5)

- Feedwater stored in Reservoir :

1. Cut the feedwater tubing (C) to the desired length (< 3 m).
2. Remove the protective plug (B) from the INLET connection (D) by pressing on collet (A) and pulling plug.
3. Connect feedwater tubing (C) to the INLET connection by inserting the tube firmly into the fitting. Verify that the connection is good by pulling back several times on the tube.
4. Connect the tube to the valve on the outlet from the reservoir. This valve should already be fitted with an adaptor to accept 8 mm diameter tubing.

- Feedwater delivered from a pressurized line in a distribution loop.

1. If the pressure is > 1 bar, then the pressure should be reduced. Use Pressure Regulator (Cat. N° ZFMQ000PR).
2. Connect feed tubing (C) to the outlet of the pressure regulator.

- **Installation of purification pack - Q-PAK
(Figure 6)**

1. Swing the pack locking handle to the left by pressing on the locking tab and pulling down on the handle. Remove the protective plugs (A).
2. Remove the protective caps (B) on the Q-PAK.
3. Wet the two O-rings (C) of the pack with a few drops of water.
4. Position the Q-PAK in its cavity. Hold the MILLI-Q 185 PLUS (D) firmly and push home the Q-PAK (E), if necessary use a slight side to side movement.
5. Replace and lock the pack retaining arm (F) in position.

3 OPERATION

• Standard modes of operation

The MILLI-Q 185 PLUS displays several modes of operation. Certain can be activated via the keyboard, others are preprogrammed and automatic.

OPERATING MODE DISPLAY	ACTIVATION	STATUS OF MILLI-Q 185 PLUS
STANDBY	Press OPERATE/STANDBY button, when system is in OPERATE mode	The system is in the STANDBY mode and is not producing purified water.
TEMPERATURE : 25°C	Press button TEMPERATURE when the system is in the STANDBY mode.	Display of the reference TEMPERATURE 25°C (Temperature compensated measurement at 25°C).
TEST : .. ; . MΩ.cm	Press button OPERATE/STANDBY when system is in STANDBY mode	Display of automatic self test result, to verify system calibration (15 ± 1 MΩ.cm displayed for 4 secs.).
PRODUCT: .. ; .MΩ.cm	Display appears automatically 4 secs after having pressed the OPERATE/STANDBY button, when the system has been in STANDBY mode. This display comes after the Auto-test display.	Display of the product water resistivity.
TEMPERATURE : ..°C	Press the button TEMPERATURE while the system is in the OPERATE or RECIRCULATION mode.	Display of the product water temperature at the Q-PAK outlet.
RECIRCULATION	Automatic function when system is in STANDBY mode.	Periodic function of the system during the STANDBY mode. The system will recirculate water automatically for 15 minutes every hour, in order to maintain water quality during periods of non-use. The 15 mn RECIRCULATION period begins 45 mn after the system is placed in STANDBY mode, and thereafter once every hour while the unit is in STANDBY. In addition, the RECIRCULATION mode is also activated each time the unit is switched ON.

• Additional Displays

DISPLAY	STATUS OF MILLI-Q 185 PLUS
PRODUCT : XX MΩ.cm	This display indicates a resistivity value which is out of scale. See TROUBLESHOOTING chapter.
TEST INCORRECT	Displayed when the Auto-test result displayed is < 14 MΩ.cm or > 16 MΩ.cm See TROUBLESHOOTING chapter.
BATTERY	The battery on the P.C. board is discharged. The data stored when the system was switched OFF has now been lost. See TROUBLESHOOTING chapter.
EXCHANGE PACK	- Flashes when the Q-PAK should be changed. See MAINTENANCE chapter. - Continously illuminated, when the Q-PAK is not in place, or incorrectly fitted. See TROUBLESHOOTING chapter.
U.V.	- Flashes for 30 seconds when the U.V. lamp is first switched ON (heating up of lamp). - Continously illuminated when the U.V. lamp is lit and its operating efficiency is at an optimum level.

• Use of Keyboard

The keyboard allows the user to control the operating mode of the system, and to display a number of preset parameters governing the correct function of the system.

DESIGNATION OF BUTTON	FUNCTION
OPERATE/STANDBY	Pressing on this button will successively activate the OPERATE and STANDBY modes.
TEMPERATURE	- During STANDBY mode : pressing this button during STANDBY mode will cause the reference temperature 25 °C to be displayed (temperature compensated measurement at 25 °C). When the button is subsequently released the temperature display will disappear after 3 seconds. - During RECIRCULATION mode or OPERATE mode : pressing this button during these modes, the current product water temperature in the system is displayed when the button is subsequently released the temperature display will disappear after 3 seconds.

• Electrical connections (at the back of system)

REMOTE STANDBY socket. This inlet allows the MILLI-Q 185 PLUS to be controlled from an external remote unit. (REMOTE STANDBY/OPERATE modes).

Contact Millipore Technical Services for further information concerning the use of this function.

• Startup (Figure 7)

1. Verify that the Q-PAK is correctly installed and that all connections have been correctly made.
2. Open the feed water isolating valve.
3. Put the point of use 3 way valve (A) (RECIRCULATION/PRODUCTION) in the PRODUCTION position.
4. Position a container (B) under the outlet fitting (C) for the MILLIPAK 40 filter unit.
5. Switch the system ON by pressing the main ON/OFF switch at the rear (see figure 1(P)).
The display should now show RECIRCULATION.
If BATTERY appears, switch the system OFF for two seconds, then back ON.
6. If the system is fed from a pressurized line via a pressure regulator, adjust the exit pressure from the regulator between 0.5 and 1 bar when the MILLI-Q 185 PLUS is in production (3 way valve (A) open in the PRODUCTION position and water flowing out of fitting (C)).
7. Leave the system running allowing water to flow out of fitting (C). This will allow all the air to be purged from the system and the Q-PAK.

Control the filling of the recipient by opening and closing the three way valve (A) a few times.

8. After having drawn 5 - 7 litres of water via fitting (C), turn the 3 way valve (A) to the RECIRCULATION position.
Press the OPERATE/STANDBY button twice to put the system into STANDBY mode. If possible, leave the system like this overnight, to allow the Q-PAK contents to fully rehydrate.
9. The following morning, press the OPERATE/STANDBY button, TEST : . . . , .MΩ.cm should be displayed for 4 seconds, then 2 seconds later the UV indicator light should flash for 30 seconds, while the lamp heats up. At the end of this 30 second period the indicator light should be constantly illuminated. Draw off 2-3 litres of water.
Put the system back into STANDBY mode by pressing once on the OPERATE/STANDBY button.

• Installation of final filter-MILLIPAK 40 (0.22 µm) (Figure 8)

When the Q-PAK has been completely purged of air, and allowed to rehydrate overnight you can proceed to fit the MILLIPAK final filter.

1. Cover the thread (A) of the MILLIPAK filter unit with 2-3 turns of Teflon tape, in order to assure a watertight seal.
2. Screw the MILLIPAK filter unit into the female fitting (C) at the point of use on the MILLI-Q 185 PLUS. Make sure that the vent (D) is accessible.
3. Press the OPERATE/STANDBY button to put the system in OPERATE mode. Turn the 3 way valve to the PRODUCTION position.

4. Purge the air contained in the MILLIPAK 40, with water, by slightly unscrewing the vent cap (E)...
5. When water exits the vent, and the air pocket has disappeared, retighten the vent cap.
6. Place the 3 way valve in the RECIRCULATION position and put the system in the STANDBY mode.
7. Insert the MILLIPAK filter into the clip (F) of the point of use cover (G).
8. Install the Millipak filter and the point of use cover on the MILLI-Q 185 PLUS by lightly pressing the tabs (H) situated on each side of the cover and inserting these tabs in the two grooves (J) intended for this purpose.
The MILLI-Q 185 PLUS is now ready for use.

4. To stop drawing water, turn the 3 way valve back to the RECIRCULATION position. Press the OPERATE/STANDBY button to put the system in the STANDBY mode.

NOTE : When you do not require water, leave the MILLI-Q 185 PLUS in the STANDBY mode.

In order to maintain water quality during non use periods the system will automatically recirculate water for 15 minutes every hour during STANDBY mode. (The display will show RECIRCULATION during this period).

During STANDBY mode make sure that the 3 way valve is in the RECIRCULATION position, to avoid water flowing out of the point of use during the 15 min RECIRCULATION.

• Withdrawing water

1. Press the button OPERATE/STANDBY (figure 1 (I)) to put the system in the OPERATE mode.
2. After 6 seconds the U.V. indicator light (figure 1 (G)) should flash. Wait a further 30 seconds until the indicator light is continuously illuminated and the resistivity of the water has reached 18.2 MΩ.cm, before drawing water. This delay allows the lamp to reach its optimum temperature.
3. Withdraw water now by turning the 3 way valve (figure 1 (D)) to the PRODUCTION position.

4 MAINTENANCE

• Maintenance schedule

FREQUENCY	ACTION
When the indicator light EXCHANGE PAK flashes, or when the resistivity no longer reaches 18.2 MΩ.cm after a few minutes of recirculation.	Replace the purification pack Q-PAK, and the final filter MILLIPAK 40.
Once a year, or when the U.V. indicator light no longer lights up.	Replace the U.V. lamp.
When the flow rate through the final filter becomes unacceptably low.	Replace the MILLIPAK 40 final filter. If this has only been in place a short time check the quality of the feedwater to the MILLI-Q 185 PLUS unit.

• Replacement of Q-PAK purification element

1. Put the MILLI-Q 185 PLUS system in STANDBY mode.

DO NOT SWITCH OFF THE SYSTEM.

2. Unlock the Q-PAK retaining arm (see figure 1(F))

3. Pull out the used Q-PAK

4. Install the replacement Q-PAK as described in INSTALLATION chapter.

NOTE : It is preferable to carry out this replacement in the afternoon or evening in order to allow the Q-PAK to rehydrate, in the system, overnight.

• Replacement of final filter MILLIPAK 40 (Figure 8)

1. Put the MILLI-Q 185 PLUS on STANDBY mode.

2. Press on the twin tabs (H) of the point of use cover (G) to remove the cover and the MILLIPAK.

3. Pull on the cover to remove it from MILLIPAK.

4. Unscrew used MILLIPAK from its fitting.

5. Install new MILLIPAK as described in INSTALLATION chapter.

• Replacement of the U.V. lamp (Figure 9)

The rear cabinet of the unit should be removed. When the MILLI-Q 185 PLUS is wall mounted, the system should be lifted down from its support bracket.

1. Switch OFF the system and remove the power cable (figure 1 (N) (P))

2. Unscrew the knurled knob (figure 1 (L)) and also the pin for the wall mounting bracket (figure 1 (K)) (if wall mounted).
3. Remove the rear cover by sliding it backwards.
4. Disconnect connector (A) from the U.V. lamp unit by pulling up.
5. Unscrew locking ring (B) withdraw the quartz tube (D) the used U.V. lamp and the O-ring (E) as shown.
6. Put in place the new quartz tube containing the new U.V. lamp.
USE GLOVES WHEN HANDLING THESE NEW ELEMENTS.
The tube should rest on a spring which will automatically centre it in the correct position.
7. Replace the O-ring (E), place the pressure ring (C) in position and screw down the locking ring (B) to tighten the whole assembly.
8. Connect the U.V. lamp (connector (A)).
9. Replace the rear cabinet, plug in the power cable, and switch the system ON.
Verify the correct functioning of the U.V. lamp : U.V. indicator light should be illuminated during OPERATE mode.

5 TROUBLESHOOTING

• Troubleshooting table

If the MILLI-Q 185 PLUS unit is not functioning correctly please look initially at the table below.
This table is an initial aid in helping you to diagnose the reason for the incorrect functioning of the system.
If this is not sufficient please contact your local Millipore Service Department. (List of addresses at the back of this manual).

PROBLEM	CAUSE	REMEDIAL ACTION
No power to MILLI-Q 185 PLUS The power light is OFF.	- Main socket is not receiving power. - Power cord not properly connected. - Main power fuse blown. - P.C. board fuse blown.	- Check mains. - Check connections. - Replace main fuse (appendix 1). - Replace P.C. board fuse (appendix 2).
EXCHANGE PACK light flashing.	Q-PAK needs to be changed.	Replace Q-PAK. See MAINTENANCE chapter.
EXCHANGE PACK light constantly illuminated.	Q-PAK is not in place, or has been incorrectly installed.	Replace Q-PAK. See MAINTENANCE chapter.
The U.V. indicator light does not light up when the system is in production or RECIRCULATION mode.	- The U.V. lamp fuse is blown. - The U.V. lamp is defective. - TEMPERATURE > 35 °C is displayed.	- Replace the fuse on the P.C. board, for the U.V. lamp (see appendix 3). - Replace the U.V lamp see MAINTENANCE chapter. - Internal water temperature is too high. See section TEMPERATURE > 35 °C indication later in this table.
BATTERY is displayed on the screen.	Battery discharged.	Leave the system under power for 48 hours to recharge the battery, switch off the system for 2 seconds to erase the display BATTERY at the end of the recharge.

PROBLEM	CAUSE	REMEDIAL ACTION
TEST INCORRECT is displayed after the autotest.	- Poor calibration of temperature circuit. - Poor calibration of resistivity circuit.	With system in STANDBY mode press button TEMPERATURE. TEMPERATURE 25 °C should be displayed. Call MILLIPORE Technical Service if the value is different. If the temperature is correctly calibrated, pressing OPERATE/STANDBY button should cause TEST : 15.0 MΩ.cm (±1 MΩ.cm) to be displayed. Call MILLIPORE Technical Service if this value is out of specification.
PRODUCT : XX MΩ.cm displayed	- Presence of air in the measurement circuit. - Fault in calibration of measurement circuit (TEST INCORRECT).	- Purge the air from Q-PAK and measurement circuit by opening the point of use valve. - Calibrate measurement circuit. (See above).
TEMPERATURE > 35 °C is displayed during PRODUCTION or RECIRCULATION mode.	- The water temperature in the system has reached 36 °C, and needs to come down to 35 °C to relight the U.V. lamp This is caused by leaving the system on OPERATE mode for too long, without drawing water. - Feed water temperature is too high.	- Put the system in STANDBY mode when not in use. - Withdraw water from the system to help to decrease the temperature, to allow the U.V. lamp to relight. - If it is not possible to reduce the feedwater temperature, then remove jumper "E" on the P.C. Board to cancel this function. Contact MILLIPORE Service Department.
TEMPERATURE < 0 °C or TEMPERATURE > + °C is displayed during OPERATE mode.	- Feedwater temperature is > 40 °C. - The thermistor, or its connection to the resistivity cell, is defective.	- Check feed water temperature and decrease it if it is elevated. - Call MILLIPORE Service department.

PROBLEM	CAUSE	REMEDIAL ACTION
The resistivity of water in the system does not reach 18.2 MΩ.cm after several minutes of recirculation. No water produced in OPERATE mode. Low flow rate of water in OPERATE mode.	The Q-PAK is exhausted. - Isolating valve on feedwater is closed. - Faulty inlet solenoid valve. - Pump running dry. - Lack of feedwater. - Air pocket build up in MILLIPAK 40 final filter. - MILLIPAK 40 final filter plugged.	Replace the Q-PAK. See MAINTENANCE chapter. - Open the feed isolating valve. - Call MILLIPORE Service Department. - The pump can pull in water only when it is primed. Introduce water into the system by gravity feed, or by pressure. - Verify the feedwater line to the system. Verify correct pressure, tubing diameter, tubing length, and flow. - Purge air from MILLIPAK 40 final filter (see OPERATION chapter, START-UP section) - Replace MILLIPAK 40 final filter. See MAINTENANCE chapter.

6 APPENDIX

• Appendix 1 : Main fuse replacement (Figure 10)

1. Switch OFF power to the MILLI-Q 185 PLUS.
2. Unplug and disconnect the main power cord to the system.
3. Remove the fuse holder (A) (figure 1 (O)).
4. Remove the faulty fuse (B) and replace by new fuse (C).

NOTE : When changing a faulty fuse, we suggest a new fuse be kept as a spare in the secondary position.

5. Replace the fuse holder in position.

• Appendix 2 : Main P.C. Board, fuse replacement (Figure 11)

1. Switch OFF power to the MILLI-Q 185 PLUS system and unplug the power cable.
2. Unscrew the knurled knob (figure 1 (L)) to remove rear cabinet.
3. Change the fuse (C) of the Main P.C. Board (A).

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.....

• Technical assistance

For further information on Millipore products, please contact your nearest Millipore subsidiary.

AUSTRALIA

Millipore Australia Pty. Limited
A.C.N. 001 239 818
Private Bag No. 18
87-89 Mars Road
Lane Cove, N.S.W.2066
Tel. Toll Free (008) 222-111
In Sydney Area (02) 428-7333
FAX (02) 427-0611

AUSTRIA, CENTRAL EUROPE, C.I.S., AFRICA, MIDDLE EAST and GULF

Millipore Ges.m.b.H.
Hietzinger Hauptstraße 145
A-1130 Wien, Austria
Tel. (43) 1 877 8926
Telex 131464 miliv a
FAX (1) 877 1654

BALTIC REPUBLICS

Millipore Oy
Ruukinkuja 1
SF-02320 Espoo 32, Finland
Tel. (90) 801 90 77
FAX (90) 801 46 40

BELGIUM and LUXEMBOURG

Millipore S.A.-N.V.
Rue de la Fusée
60/Raketstraat 60
B-1130 Bruxelles-Brussel
Tel. (02) 242.17.40
FAX (02) 242.27.85

BRAZIL

Millipore Indústria e Comércio
Ltda.
Rua Prof. Campos de Oliveira,
430
CEP 04675
São Paulo-SP
Tel. (011) 548-7011
Telex 11 57751 MILY BR
FAX (011) 548-7923

CANADA

Millipore (Canada) Ltd.
3688 Nashua Drive
Mississauga, Ontario/L4V 1M5
Tel. Toll Free 1-800-268-4881
In Toronto Area (416) 678-2161
FAX (416) 678-0882

CHINA, People's Republic of

Millipore China Ltd.
Suite 1101-1106 and 16
Asia Pacific Bldg.
No. 8 Ya Bao Road
Chao Yang District
Beijing
Postal Code 100020
Tel. (86-1) 5135114/5135116-
7/5125992/5126123
Telex 22685 NMLBJ
FAX (86-1) 5126113

Millipore China Ltd.
Room 615, China Hotel Office
Tower
Linhua Lu
Guangzhou
Postal Code 510015
Tel. (86) 20-686217

Millipore China Ltd.
Suite 2805, Union Building
100 Yan An Road (E)
Shanghai 200002
Tel. (86) 21-3203856
FAX (86) 21-3200236

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• Warranty

Millipore Corporation (Millipore) warrants the products manufactured by it against defects in materials and workmanship when used in accordance with the applicable instructions for a period of one year from the date of shipment of the products. MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

The warranty provided herein and the data, specification and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period, provided the customer notifies Millipore promptly of any such defect. The exclusive remedy provided herein shall not be deemed to have failed of its essential purpose so long as Millipore is willing to repair or replace any nonconforming Millipore product or part. Millipore shall not be liable for consequential damages resulting from economic loss or property damages sustained by a customer from the use of its products.

However, in some states the purchaser may have rights under state law in addition to those provided by this warranty.

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Millipore (U.K.) Ltd.
U.S.A.